


# Hypotube



Ø - tw - x - l / Hypotube with stiletto (lasered)



Ø - tw - x - l / Hypotube with skive



Ø - tw - x - l / Hypotube (laser cutting)

## Explanations

- Ø diameter of the tube  
(standard 0.3 - 2.0; others on request)
- tw wall thickness in mm (from 0.05)
- x coatings
  - U uncoated
  - P PTFE
  - H hydrophil
- l lengths (individually on request)

**Material:** stainless steel 1.4301

**Coating:** individually on request

**Tolerances:** depends on the length  
 +/- 0.2 mm bis +/- 1.0 mm

## The different possibilities of processing are:

- sandblasted
- slotted
- welded (laser or microplasma)
- ground
- laser cutting
- others on request

## Customer requirements

Besides our standards, we would be pleased to be your manufacturer when it comes to individual specifications and wishes for Hypotubes.