


## Niro-Stonecatcher Basket



NiRo 4 (6) - DØ - KØ

### Construction

On the new Niro-Stone catcher are all transitions soldered, which guarantees a durable and firm connection, which makes the Basket lithotripsy capable.

### Explanations

DØ- diameter of wire  
 KØ- Basket Diameter  
 (20, 25, 30, 35 mm)

### Technical Data

#### Basket

4 (Standard) or 6 stainless steel wire ropes summarized in a twisted stainless steel head.

#### Rope

Consisting of a stainless steel wire rope. After requirement a stainless steeldrawbar can be mounted.

### Customer requirements

Besides our standards, we would be pleased to be your manufacturer when it comes to individual specifications and wishes for Niro-Stonecatcher Baskets.