



UROLOGY



endox®
Feinwerktechnik GmbH



EXTRACTOR BASKETS (NiTi)



TIPLESS


options	120 cm	90 cm
length	120 cm	90 cm
basket Ø	12 mm	14 mm
catheter Ø	1.9 Fr	3.0 Fr
cathetermaterial	Polyimide	PEEK
handle	3-finger	
features	360° free rotatable completely deattachable single wire technology - braided basket without connecting points	

WITH TIP

options	90 cm	120 cm
length	90 cm	120 cm
basket Ø	10 mm	bis 15 mm
catheter Ø	1.9 Fr	3.0 Fr
cathetermaterial	Polyimide	PEEK
forms	diamond spherical shape double helical helical	




360° free rotatable




complete detachable




single wire technology



spherical shape



double helical



helical


GUIDE WIRES (NiTi)

PU-EXTRUDED/HYDROPHILIC

specification	
Ø	0,035 inch
PU	PU-extruded
L	150 cm
J	60° (option)
F	flexible on one side

PTFE WITH PU-TIP

specification	
Ø	0,035 inch
PU	PU-tip
P	PTFE
L	150 cm
J	60° (option)
F	flexible on one side
FF	flexible on both sides



Customer requirements

In addition to our standards, we manufacture according to your specifications and wishes.