

HYPOTUBE



Ø - tw - x - l / Hypotube with stylet
(lasered)



Ø - tw - x - l / Hypotube with skive



Ø - tw - x - l / Hypotube
(laser cutting)

The different possibilities of processing are:

- sandblasting
- slitting
- welding (laser or microplasma)
- grinding
- laser cutting
- others on request

Customer requirements

In addition to our standards we also manufacture hypotubes according to your specifications and wishes.

Explanations

Ø	diameter of the tube (standard 0.3 - 2.0; others on request)
tw	wall thickness in mm (from 0.05)
x	coatings U uncoated P PTFE H hydrophilic
l	lengths in mm (individually up to 5.000)
Material	stainless steel 1.4301
Coatings	individually on request
Tolerances	depends on the length +/- 0.2 mm bis +/- 1.0 mm

