

# STAINLESS STEEL STONECATCHER



NiRo 4 (6) - DØ - KØ

## Construction

All transitions of the stainless steel stone catcher are soldered, which ensures a durable and strong connection that makes the basket suitable for lithotripsy.

## Technical Data

*Basket* 4 (standard) or 6 stainless steel wire ropes summarized in a twisted stainless steel tip.

*Rope* Made of a stainless steel strand rope. If necessary, a stainless steel pull rope can be attached.

## Explanations

DØ diameter of wire  
KØ diameter of basket  
(20, 25, 30, 35 mm)

## Customer requirements

In addition to our standards we also manufacture stainless steel stonecatchers according to your specifications and wishes.

