

WIRE / STRAND SECTIONS



Ø - x - l
Ø - l - k



In addition to our standards we also manufacture wire- / strand- sections according to your specifications and wishes.

Explanations (wire)

Ø	diameter in mm (standard 0.1 - 1.0 mm; other on request)
x	coatings U uncoated P PTFE H hydrophilic
l	lengths in mm (individually on request up to 5.000)
Material	stainless steel 1.4310 or Nitinol
Coating	individual
Taper	individual
Tolerances	depends on the length +/- 1 mm up to +/- 3 mm

Explanations (strand)

Ø	nominal diameter in mm (standard 0.09 bis 2.0)
l	lengths in mm (bis 5.000)
k	type of construction (1 x 7, 1 x 19, 7 x 7, und 7 x 19)
Material	stainless steel 1.4401 or 1.4310
End finish	rounded, balls are individually possible
Tolerances	depends on length +/- 1 mm up to +/- 3 mm

